

3SHAPE COMMUNICATE

Connectie Elysee Dental




3shape 

 **ELYSEEDENTAL**
Modern Dental Group

CONNECTIE

- 1 www.3shapecommunicate.com
- 2 Klik op Connections → Add Connection



- 3 scandata@elysee-dental.com → Show me Labs → Find → Connect

A screenshot of the 'Connect to a Practice or Lab' form. The title is 'Connect to a Practice or Lab'. Below the title, there is a text input field with the email 'scandata@elysee-dental.com' and a 'Find' button. Below the input field, there are two radio buttons: 'Show me Labs' (checked) and 'Show me Practices'. A search result for 'Elysee Dental BE' is displayed, including the email 'scandata@elysee-dental.com', location 'Leuven, BE', address 'Minderbroederstraat 12Q', phone number '0032 16 27 00 40', and a 'TRIOS Ready' logo. Below the search result, there is a text area for a message to the practice, containing the text 'We would like to connect with you on 3Shape Communicate.' and a 'Connect' button.

- 4 Uw aanvraag wordt zo snel mogelijk in behandeling genomen